

Technical Data Sheet

Product **Anti-GlcB (M. tuberculosis) Purified**

Cat. Number/Size **11-192-C100 0.1 mg**

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	GlcB (M. tuberculosis)
Clone	PAb (192)
Format	Purified
Reactivity	Mycobacterium tuberculosis
Application	WB, ELISA
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody to GlcB reacts with Mycobacterium tuberculosis protein GlcB (malate synthase, Rv1837c).
Other names	Rv1837c
Immunogen	Recombinant GlcB protein (produced in E. coli)
Entrez Gene ID	885713
Gene name	glcB
NCBI Full Gene Name	malate synthase
UniProt ID	P9WK17

Concentration	1 mg/ml
Preparation	Purified from rabbit serum by affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.