



	Technical Data Sheet	
Product	Anti-Hu CD17 Purified	
Cat. Number/Size	11-666-C025	0.025 mg
	11-666-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD17	
Clone	MEM-68	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2 µg/ml.	
Isotype	Mouse IgM	
Specificity	The mouse monoclonal antibody MEM-68 recognizes CD17, an extracellular membrane lipid moiety lactosylceramide expressed on granulocytes, monocytes and platelets.	
Other names	Lactosylceramide	
Immunogen	Pre-B cell line NALM-6	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.