



	Technical Data Sheet	
Product	Anti-HLA-G Purified	
Cat. Number/Size	11-432-C025	0.025 mg
	11-432-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HLA-G	
Clone	2A12	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), WB, IHC-P, IHC-F, ELISA	
Application details	Flow cytometry: Excellent for intracellular staining; recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody 2A12 was generated to a peptide corresponding to C-intron 4-encoded sequence. This antibody does not crossreact with the full-length HLA-G1 isoform and thus allows to distinguish between secreted HLA-G5 and HLA-G6 isoforms from shedded HLA-G1.	
Immunogen	C-terminal amino acid sequence (22-mer) of soluble HLA-G5 and HLA-G6 proteins coupled to ovalbumin.	
Entrez Gene ID	3135	
Gene name	HLA-G	
NCBI Full Gene Name	major histocompatibility complex, class I, G	
UniProt ID	P17693	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.