



	Technical Data Sheet	
Product	Anti-Ms CD49e Purified	
Cat. Number/Size	11-743-C025	0.025 mg
	11-743-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD49e	
Clone	5H10-27 (MFR5)	
Format	Purified	
Reactivity	Mouse	
Application	ICC, FC (QC tested), IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Rat IgG2a kappa	
Specificity	The rat monoclonal antibody 5H10-27 (MFR5) recognizes an extracellular epitope of murine CD49e, a 135 kDa protein serving as VLA-5 alpha chain, expressed on thymocytes, activated T cells, splenic B cells, and mast cells.	
Other names	ITGA5, VLA-5 alpha, FNRA	
Immunogen	C57BL/6 x A/J)F1 murine mast cell line	
Entrez Gene ID	16402	
Gene name	Itga5	
NCBI Full Gene Name	integrin alpha 5	
UniProt ID	P11688	
Concentration	1 mg/ml	
Preparation	Purified by protein-G affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.