

Technical Data Sheet

Product

Anti-Hu CD2 FITC

Cat. Number/Size

1F-492-T025

25 tests

1F-492-T100

100 tests

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD2
Clone	LT2
Format	FITC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG2b
Specificity	The antibody LT2 reacts with an extracellular epitope of CD2, a 50 kDa glycoprotein present on the human peripheral blood T lymphocytes and NK cells; also expressed by all thymocytes.
Other names	T11, LFA-2, SRBC
Workshop	HLDA VI: WS Code T 6T-008
Immunogen	Normal human blood lymphocytes.
Entrez Gene ID	914
Gene name	CD2
NCBI Full Gene Name	CD2 molecule
UniProt ID	P06729

Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.