



Technical Data Sheet

Product	Anti-Hu CD11b Purified Azide Free	
Cat. Number/Size	10-502-C025	0.025 mg
	10-502-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		

Antigen	CD11b
Clone	MEM-170
Format	Purified Azide Free
Reactivity	Human
Application	FC (QC tested), FUNC
Application details	Functional application: The antibody MEM-170 blocks (not completely) Mac-1-mediated neutrophil binding to fibrinogen. Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody MEM-170 recognizes an extracellular epitope of CD11b antigen (Mac-1alpha), a 165-170 kDa transmembrane protein type 1 mainly expressed on monocytes, granulocytes and NK-cells. The antibody blocks (not completely) Mac-1-mediated neutrophil binding to fibrinogen.
Other names	Mac-1, Integrin alpha M, ITGAM, CR3A, MO1A, MAC1A
Workshop	HLDA VI: WS Code A054
Immunogen	Human granulocytes
Entrez Gene ID	3684
Gene name	ITGAM
NCBI Full Gene Name	integrin subunit alpha M
UniProt ID	P11215
Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.
Images and References	www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.