



Technical Data Sheet

Product **Anti-Hu Ig Kappa Light Chain Purified**

Cat. Number/Size **11-674-C025** **0.025 mg**
11-674-C100 **0.1 mg**

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	Ig Kappa Light Chain
Clone	TB28-2
Format	Purified
Reactivity	Human
Application	FC (QC tested), WB
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Extracellular and intracellular staining.
Isotype	Mouse IgG1 kappa
Specificity	The antibody TB28-2 recognizes human Ig kappa light chains of both secreted and cell surface immunoglobulin. It detects also free kappa light chains.
Other names	Immunoglobulin kappa, Igk
Immunogen	Human IgG-kappa myeloma protein
Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.
Images and References	www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.