



	Technical Data Sheet	
Product	Anti-Hu CD85j Purified	
Cat. Number/Size	11-740-C025	0.025 mg
	11-740-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD85j	
Clone	GHI/75	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG2b kappa	
Specificity	The mouse monoclonal antibody GHI/75 recognizes an extracellular epitope of CD85j / ILT2, an 110-120 kDa membrane glycoprotein expressed strongly on plasma cells, moderately on circulating B cells, and weakly on monocytes. It is also expressed on T cell and NK cell subsets (variable, individual).	
Other names	ILT2, LIR1, MIR7, PIRB, ILT-2	
Immunogen	Hairy cell leukaemia cells	
Entrez Gene ID	10859	
Gene name	LILRB1	
NCBI Full Gene Name	leukocyte immunoglobulin like receptor B1	
UniProt ID	Q8NHL6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.