



	Technical Data Sheet	
Product	Anti-Hu CD41a Purified	
Cat. Number/Size	11-417-C025	0.025 mg
	11-417-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD41a	
Clone	HIP8	
Format	Purified	
Reactivity	Non-human primates, Human	
Application	FC (QC tested), IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1	
Specificity	The antibody HIP8 reacts with an extracellular epitope on alpha (a) subunit of CD41 (heavy chain; 120 kDa). CD41 is mainly expressed on platelets and megakaryocytes. The antibody HIP8 blocks platelet aggregation and completely inhibits ADP-, epinephrine-, and collagen-induced platelet activation, and partially inhibits ristocetin- and thrombin-induced platelet activation. HIP8 is useful in the morphological and physiological studies o platelets and megakaryocytes.	
Other names	GT, GTA, GP2B, HPA3, GPIIb	
Workshop	HLDA IV: WS Code P 38	
Entrez Gene ID	3674	
Gene name	ITGA2B	
NCBI Full Gene Name	integrin subunit alpha 2b	
UniProt ID	P08514	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.