



	Technical Data Sheet	
Product	Annexin V Biotin	
Cat. Number/Size	EXB0031	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Format	Biotin	
Application	FC (QC tested)	
Description	Annexin V Biotin is intended for detection of early apoptotic cells by flow cytometry.	
Contents	Annexin V Biotin - 1 vial (0.2 ml).	
Concentration	0.5 mg/ml	
Preparation	Purified recombinant protein is conjugated with biotin LC-NHS under optimum conditions. Unconjugated protein and free biotin are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.