

Technical Data Sheet

Product	Anti-Ms CD11a Purified Low Endotoxin	
Cat. Number/Size	12-622-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD11a
Clone	M17/4
Format	Purified Low Endotoxin
Reactivity	Mouse
Application	FC (QC tested), IP, IHC(F), FUNC
Application details	Functional application: Blocking. Flow cytometry: Recommended dilution: 1 µg/ml. Immunohistochemistry (frozen sections): Positive tissue: murine spleen or thymus, acetone fixation.
Isotype	Rat IgG2a
Specificity	The rat monoclonal antibody M17/4 reacts with an extracellular epitope of CD11a (alpha-subunit of murine LFA-1), a 180 kDa type I transmembrane glycoprotein expressed on B and T lymphocytes, monocytes, macrophages, neutrophils, basophils and eosinophils.
Other names	LFA-1, LFA1A, ITGAL
Immunogen	C57BL/6 mouse splenic secondary cytotoxic T lymphocytes
Entrez Gene ID	16408
Gene name	Itgal
NCBI Full Gene Name	integrin subunit alpha L
UniProt ID	E9Q5M7

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.