

Technical Data Sheet

Product

**Anti-ARAP1 Purified**

Cat. Number/Size

**11-516-C025**

**0.025 mg**

**11-516-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	ARAP1
Clone	ARAP1-2
Format	Purified
Reactivity	Human
Application	IP, WB, ICC
Application details	Western blotting: Recommended dilution: 2 µg/ml
Isotype	Mouse IgG1
Specificity	The antibody ARAP1-2 reacts with C-terminal part of ARAP1 (intracellular epitope), a 160 kDa adaptor protein.
Other names	CED2, CENTD2
Immunogen	His6-tagged protein corresponding to amino acids 1190-1450 of human ARAP1
Entrez Gene ID	116985
Gene name	ARAP1
NCBI Full Gene Name	ArfGAP with RhoGAP domain, ankyrin repeat and PH d
UniProt ID	Q96P48

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.