

Technical Data Sheet

Product	<b>Anti-Hu CD88 Purified</b>	
Cat. Number/Size	<b>11-736-C025</b>	<b>0.025 mg</b>
	<b>11-736-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD88
Clone	S5/1
Format	Purified
Reactivity	Rabbit, Ferret, Cow, Human
Application	FC (QC tested), IP, WB, IHC(P)
Application details	Flow cytometry: Recommended dilution: 3-12 µg/ml.
Isotype	Mouse IgG2a kappa
Specificity	The mouse monoclonal antibody S5/1 recognizes an extracellular epitope of CD88 protein, a 43 kDa receptor of C5a component of the complement cascade.
Other names	C5A, C5R1, C5AR1
Workshop	HCDM listed
Immunogen	Recombinant N-terminal peptide (Asp15-Asp27) of human C5aR
Entrez Gene ID	728
Gene name	C5AR1
NCBI Full Gene Name	complement C5a receptor 1
UniProt ID	P21730

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.