

Technical Data Sheet

Product	<b>Anti-Hu LAR Purified</b>	
Cat. Number/Size	<b>11-830-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	LAR
Clone	W7C6
Format	Purified
Reactivity	Non-human primates, Human
Application	FC (QC tested)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody W7C6 recognizes an extracellular epitope of protein tyrosine phosphatase LAR, a marker of mesenchymal stem cells.
Other names	Leucocyte common antigen related molecule, PTPRF, W7C6
Immunogen	WERI-RB-1 retinoblastoma cells
Entrez Gene ID	5792
Gene name	PTPRF
NCBI Full Gene Name	protein tyrosine phosphatase receptor type F
UniProt ID	G1UI20

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.