

Technical Data Sheet

Product	<b>Anti-Cytokeratin 10 Purified</b>	
Cat. Number/Size	<b>11-111-C025</b>	<b>0.025 mg</b>
	<b>11-111-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	Cytokeratin 10
Clone	VIK-10
Format	Purified
Reactivity	Human
Application	IHC(P), IP, ICC
Isotype	Mouse IgG1
Specificity	The antibody VIK-10 reacts with Cytokeratin 10 (56.5 kDa; intracellular antigen). Cytokeratins are members of intermediate filaments subfamily represented in epithelial tissues.
Other names	BIE, EHK, K10, KPP, BCIE, CK10, CK-10, CYK10, KRT10
Immunogen	Cytoskeleton preparation extracted from human epidermis by detergent/high salt extraction.
Entrez Gene ID	3858
Gene name	KRT10
NCBI Full Gene Name	keratin 10
UniProt ID	P13645

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.