



	Technical Data Sheet	
Product	Anti-Hu CD158f Purified	
Cat. Number/Size	11-774-C025	0.025 mg
	11-774-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD158f	
Clone	UP-R1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody UP-R1 recognizes an extracellular epitope on CD158f (KIR2DL5), a 60 kDa glycoprotein serving as a HLA class I ligand, and mainly expressed on a subset of NK cells and a small population of T cells. Its expression is highly polymorphic between individuals.	
Other names	KIR2DL5A	
Immunogen	Human CD158f-Ig fusion protein	
Entrez Gene ID	57292	
Gene name	KIR2DL5A	
NCBI Full Gene Name	killer cell immunoglobulin like receptor, two Ig d	
UniProt ID	Q8N109	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.