

Technical Data Sheet

Product	Anti-CD3 zeta PE
Cat. Number/Size	1P-568-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	CD3 zeta
Clone	H146-968
Format	PE
Reactivity	Mouse, Human
Application	FC-IC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Excitation laser	blue (488 nm)
Isotype	Armenian Hamster IgG2 kappa
Specificity	The Armenian hamster antibody H146-968 reacts with an intracellular epitope (amino acids 151-164) of CD3 zeta chain (CD247), which is a component of TCR/CD3 complex expressed on T cells.
Other names	CD3 zeta chain, CD247, T3Z, TCRz, IMD25
Immunogen	Synthetic peptide corresponding to amino acids 151-164 of mouse CD3 zeta.
Entrez Gene ID	919
Gene name	CD247
NCBI Full Gene Name	CD247 molecule
UniProt ID	P20963

Concentration	0.1 mg/ml
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.