



	Technical Data Sheet	
Product	Anti-Hu CD264 PE	
Cat. Number/Size	1P-519-C025	0.025 mg
	1P-519-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD264	
Clone	TRAIL-R4-01	
Format	PE	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG1	
Specificity	The antibody TRAIL-R4-01 reacts with an extracellular epitope of TRAIL-R4, a 42 kDa transmembrane protein expressed on various blood cells.	
Other names	TNFRSF10D, DcR2, TRUNDD, TRAILR4	
Workshop	HLDA IX	
Immunogen	TRAIL-R4 (aa 1-210) - hIgGhc fusion protein	
Entrez Gene ID	8793	
Gene name	TNFRSF10D	
NCBI Full Gene Name	TNF receptor superfamily member 10d	
UniProt ID	Q9UBN6	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.