



	Technical Data Sheet	
Product	Anti-HRP Purified	
Cat. Number/Size	11-262-C025	0.025 mg
	11-262-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HRP	
Clone	HP-03	
Format	Purified	
Reactivity	Horseradish, Tagged fusion proteins in all species	
Application	WB, ICC, ELISA	
Application details	Immunocytochemistry: The antibody HP-03 was successfully used for staining of formaldehyde-fixed, Triton-permeabilized cells transfected with HRP gene.	
Isotype	Mouse IgG1	
Specificity	The antibody HP-03 reacts with horseradish peroxidase. It is suitable for PAP detection system.	
Immunogen	Horseradish peroxidase (Grade I)	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.