



	Technical Data Sheet	
Product	Anti-Ms CD3 PerCP-Cy™5.5	
Cat. Number/Size	T9-578-C025	0.025 mg
	T9-578-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD3	
Clone	145-2C11	
Format	PerCP-Cy™5.5	
Reactivity	Mouse	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-3 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Armenian Hamster IgG kappa	
Specificity	The Armenian hamster monoclonal antibody 145-2C11 reacts with an extracellular epitope of murine CD3 (epsilon subunit). This antibody is commonly used as a phenotypic marker for murine T cells.	
Other names	CD3E, T3E, TCRE	
Immunogen	Mouse BM10-37 cytotoxic T lymphocytes	
Entrez Gene ID	12501	
Gene name	CD3E	
NCBI Full Gene Name	CD3 antigen, epsilon polypeptide	
UniProt ID	P22646	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with activated tandem dye of activated Peridinin-Chlorophyll Protein-Cyanine 5.5 (PerCP-Cy™5.5) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.

Cy™ and CyDye™ are registered trademarks of Cytiva.