



	Technical Data Sheet	
Product	Anti-TNF alpha Purified	
Cat. Number/Size	11-687-C025	0.025 mg
	11-687-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	TNF alpha	
Clone	MAb1	
Format	Purified	
Reactivity	Human	
Application	WB, ELISA	
Application details	ELISA: It can be used as capture antibody in combination with biotinylated antibody MAb11.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody MAb1 recognizes human 17-26 kDa cytokine TNF-alpha (tumor necrosis factor alpha).	
Other names	TNFA, DIF, TNLG1F	
Immunogen	Recombinant human TNF-alpha	
Entrez Gene ID	7124	
Gene name	TNF	
NCBI Full Gene Name	tumor necrosis factor	
UniProt ID	P01375	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.