



	Technical Data Sheet	
Product	Anti-Hu Ig Kappa Light Chain Purified	
Cat. Number/Size	11-504-C025	0.025 mg
	11-504-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Ig Kappa Light Chain	
Clone	A8B5	
Format	Purified	
Reactivity	Human	
Negative species	Rabbit, Hamster, Sheep, Guinea pig, Goat	
Application	FC (QC tested), IHC-P, IHC-F, ELISA	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml, Positive control: Daudi human cell line.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody A8B5 reacts with kappa light chains (22.5 kDa) of immunoglobulins.	
Other names	Immunoglobulin kappa, Igk	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.