

Technical Data Sheet

Product	<b>Anti-LIME1 Purified</b>	
Cat. Number/Size	<b>11-400-C025</b>	<b>0.025 mg</b>
	<b>11-400-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	LIME1
Clone	LIME-06
Format	Purified
Reactivity	Human
Application	IP, FC(IC)
Application details	Flow cytometry: Recommended dilution: 8-12 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Specificity	The antibody LIME-06 was raised against intracellular fragment corresponding to aa 141-295 of human LIME, a 30 kDa Lck-interacting transmembrane adaptor expressed by T cells.
Other names	LITA1
Immunogen	Bacterially expressed intracellular fragment corresponding to aa 141-295 of human LIME.
Entrez Gene ID	54923
Gene name	LIME1
NCBI Full Gene Name	Lck interacting transmembrane adaptor 1
UniProt ID	Q9H400

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.