



	Technical Data Sheet	
Product	Anti-Hu CD267 Purified	
Cat. Number/Size	11-162-C025	0.025 mg
	11-162-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD267	
Clone	1A1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Rat IgG2a kappa	
Specificity	The rat monoclonal antibody 1A1 recognizes an extracellular epitope of CD267 / TACI, a 32 kDa type III transmembrane protein expressed by B cells and possibly by some activated T cells.	
Other names	RYZN, TACI, CVID2, IGAD2, TNFRSF14B	
Workshop	HLDA VIII	
Immunogen	CD267-transfected RBL cells	
Entrez Gene ID	23495	
Gene name	TNFRSF13B	
NCBI Full Gene Name	TNF receptor superfamily member 13B	
UniProt ID	O14836	
Concentration	1 mg/ml	
Preparation	Purified by protein-G affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.