

Technical Data Sheet

Product	Anti-Dendra2 purified
Cat. Number/Size	11-836-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	Dendra2
Clone	PAb (836)
Format	Purified
Reactivity	Species independent
Application	WB, ICC, FC (QC tested)
Application details	Western blotting: Recommended dilution: 1-2 µg/ml. Flow cytometry: Recommended dilution: 1-4 µg/ml. Extracellular and intracellular staining.
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody PAb (836) recognizes Dendra2 tag, a green-to red photoswitchable protein derived from an octocoral (Dendronephthya).
Other names	Dendra2 tag
Immunogen	Dendra2 tag

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.