



	Technical Data Sheet	
Product	Anti-LYN Purified	
Cat. Number/Size	11-377-C025	0.025 mg
	11-377-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	LYN	
Clone	LYN-01	
Format	Purified	
Reactivity	Rat, Human, Mouse	
Application	IP, WB, ICC	
Application details	Western blotting: Recommended dilution: 1-2 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody LYN-01 reacts with Lyn (p56/p53), a non-receptor Src-family tyrosine kinase expressed in hematopoietic tissues (intracellular antigen).	
Other names	JTK8, p56Lyn, p53Lyn	
Immunogen	Bacterially expressed recombinant fragment of human Lyn (aa 8 - 238).	
Entrez Gene ID	4067	
Gene name	LYN	
NCBI Full Gene Name	LYN proto-oncogene, Src family tyrosine kinase	
UniProt ID	P07948	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.