



	Technical Data Sheet	
Product	Anti-Hu CD300a Purified	
Cat. Number/Size	11-501-C025	0.025 mg
	11-501-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD300a	
Clone	MEM-260	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), WB	
Application details	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Positive material: HPB cell line, negative material: JURKAT human cell line, RAJI human cell line; non-reducing conditions.	
Isotype	Mouse IgG1	
Specificity	The antibody MEM-260 reacts with an extracellular epitope of CD300a, a 60 kDa leukocyte transmembrane glycoprotein expressed on human granulocytes, monocytes, neutrophils, NK cells, mast cells and dendritic cells, 25% of circulating T cells and 15% of circulating B cells.	
Other names	IRC, CLM-8, IGSF12, CMRF35H	
Immunogen	HPB human acute lymphoid leukemia cell line	
Entrez Gene ID	11314	
Gene name	CD300A	
NCBI Full Gene Name	CD300a molecule	
UniProt ID	Q9UGN4	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.