

#### Technical Data Sheet

Product	<b>Anti-Hu CD73 Purified</b>	
Cat. Number/Size	<b>11-675-C025</b>	<b>0.025 mg</b>
	<b>11-675-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD73
Clone	AD2
Format	Purified
Reactivity	Human
Application	FC (QC tested), MC (CyTOF)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody AD2 recognizes CD73, a 70 kDa GPI-anchored 5'-nucleotidase expressed predominantly on the surface of T and B cell subsets, follicular dendritic cells and endothelial cells.
Other names	ecto-5'-nucleotidase, NT5E, E5NT, CALJA, NTE
Entrez Gene ID	4907
Gene name	NT5E
NCBI Full Gene Name	5'-nucleotidase ecto
UniProt ID	P21589

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.