



	Technical Data Sheet	
Product	Anti-Hu CD20 Purified	
Cat. Number/Size	11-638-C025	0.025 mg
	11-638-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD20	
Clone	2H7	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), IP, IHC-F, ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG2b kappa	
Specificity	The mouse monoclonal antibody 2H7 recognizes an extracellular epitope on CD20 (B1, Bp35), a 33-37 kDa non-glycosylated membrane receptor with four transmembrane domains, expressed on pre-B lymphocytes, resting and activated B cells (not plasma cells), follicular dendritic cells, and at low levels on peripheral blood T lymphocytes.	
Other names	B1, S7, MS4A, Bp35, CVID5, LEU-16	
Workshop	HLDA IV: WS Code B201	
Immunogen	Human tonsillar B cells	
Entrez Gene ID	931	
Gene name	CD20	
NCBI Full Gene Name	MS4A1	
UniProt ID	P11836	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.