

Technical Data Sheet

Product

**Anti-LYN Purified**

Cat. Number/Size

**11-377-C025**

**0.025 mg**

**11-377-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Antigen

LYN

Clone

LYN-01

Format

Purified

Reactivity

Rat, Human, Mouse

Application

IP, WB, ICC

Application details

Western blotting: Recommended dilution: 1-2 µg/ml.

Isotype

Mouse IgG1

Specificity

The antibody LYN-01 reacts with Lyn (p56/p53), a non-receptor Src-family tyrosine kinase expressed in hematopoietic tissues (intracellular antigen).

Other names

JTK8, p56Lyn, p53Lyn

Immunogen

Bacterially expressed recombinant fragment of human Lyn (aa 8 - 238).

Entrez Gene ID

4067

Gene name

LYN

NCBI Full Gene Name

LYN proto-oncogene, Src family tyrosine kinase

UniProt ID

P07948

Concentration

1 mg/ml

Preparation

Purified by protein-A affinity chromatography.

Formulation

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.