

Technical Data Sheet

Product	Anti-Hu IgM PerCP-Cy™5.5	
Cat. Number/Size	T9-320-C025	0.025 mg
	T9-320-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	IgM
Clone	CH2
Format	PerCP-Cy™5.5
Reactivity	Human
Application	FC (QC tested)
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The antibody CH2 reacts with Fc fragment of human IgM.
Other names	immunoglobulin M
Immunogen	Purified human IgM.

Concentration	0.1 mg/ml
Preparation	The purified antibody is conjugated with tandem dye PerCP-Cy™5.5 under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.

Cy™ and CyDye™ are registered trademarks of GE Healthcare.