

Technical Data Sheet

Product	<b>Anti-Hu STRO-1 Purified</b>	
Cat. Number/Size	<b>11-649-C025</b>	<b>0.025 mg</b>
	<b>11-649-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	STRO-1
Clone	STRO-1
Format	Purified
Reactivity	Human, Rabbit, Rat
Application	FC (QC tested), ICC
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgM lambda
Specificity	The mouse monoclonal antibody STRO-1 recognizes an extracellular epitope of the cell surface antigen STRO-1 expressed by bone marrow mesenchymal stromal cells and nucleated erythroid precursors, but not by committed hematopoietic progenitors.
Other names	STRO1
Immunogen	Human CD34 positive bone marrow cells

---

Concentration	1 mg/ml
Preparation	Purified by precipitation and chromatography
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>
-----------------------	--

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.