



	Technical Data Sheet	
Product	Anti-Hu CD272 Purified	
Cat. Number/Size	11-877-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD272	
Clone	MIH26	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2-5 µg/ml.	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody MIH26 recognizes an extracellular epitope of CD272, a transmembrane glycoprotein serving as a negative regulator of the activation in various leukocyte types.	
Other names	BTLA1	
Immunogen	cells transfected with human CD272	
Entrez Gene ID	151888	
Gene name	BTLA	
NCBI Full Gene Name	B and T lymphocyte associated	
UniProt ID	Q7Z6A9	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.