

Technical Data Sheet

Product	Anti-Benzo[a]pyrene Purified	
Cat. Number/Size	11-339-C025	0.025 mg
	11-339-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Benzo[a]pyrene	
Clone	BAP-13	
Format	Purified	
Application	ICC, ELISA	
Isotype	Mouse IgG1	
Specificity	The antibody BAP-13 is specific for Polynuclear Aromatic Hydrocarbon and inhibits DNA adduct formation (this ability was investigated in rat liver microsomes spiked with calf thymus DNA and 7,8-diol-B[a]P).	
Immunogen	Benzo[a]pyrenyl-1-butyric acid conjugated to BSA.	
Concentration	1 mg/ml	
Preparation	Purified by precipitation and chromatography	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.