

Technical Data Sheet

Product

**Anti-Hu IgM Purified**

Cat. Number/Size

**11-321-C025**

**0.025 mg**

**11-321-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Antigen

IgM

Clone

MA2

Format

Purified

Reactivity

Human

Application

ELISA

Isotype

Mouse IgG2b

Specificity

The antibody MA2 reacts with  $\mu$ -chain of human IgM, Fab - region localized epitope.

Other names

immunoglobulin M

Immunogen

IgM from human serum.

Concentration

1 mg/ml

Preparation

Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).

Formulation

Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.