



	Technical Data Sheet	
Product	Anti-RLTPR Purified	
Cat. Number/Size	11-718-C025	0.025 mg
	11-718-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	RLTPR	
Clone	EM-53	
Format	Purified	
Reactivity	Mouse, Human	
Application	FC (QC tested), WB	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody EM-53 recognizes RLTPR / CARMIL2, an intracellular protein playing a role in actin filament elongation.	
Other names	CARMIL2, LRRC16C	
Immunogen	Murine RLTPR	
Entrez Gene ID	146206	
Gene name	CARMIL2	
NCBI Full Gene Name	capping protein regulator and myosin 1 linker 2	
UniProt ID	Q6F5E8	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.