



	Technical Data Sheet	
Product	Anti-Hu CD3 Purified Low Endotoxin	
Cat. Number/Size	12-631-C100	0.1 mg
	12-631-M001	1.0 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD3	
Clone	OKT3	
Format	Purified Low Endotoxin	
Reactivity	Human	
Application	FC (QC tested), IHC-F, FUNC	
Application details	Functional application: counteracting both generation and function of effector T cells. Flow cytometry: Recommended dilution: 1 µg/ml.	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody OKT3 recognizes an extracellular epitope on CD3 antigen of the TCR/CD3 complex on mature human T cells. This antibody, also known as Orthoclone OKT3 or Muromonab-CD3, has been extensively used as a drug for therapy of acute, glucocorticoid resistant rejection of allogenic renal, heart and liver transplants. It has also been investigated for use in treating T-cell acute lymphoblastic leukemia.	
Other names	CD3E, T3E, TCRE	
Workshop	HLDA VI: WS Code T6T-R3	
Immunogen	human T cells	
Entrez Gene ID	916	
Gene name	CD3E	
NCBI Full Gene Name	CD3 antigen, epsilon polypeptide	
UniProt ID	P07766	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.