

Technical Data Sheet

Product	Anti-Hu IgG (Fab) PE-Cy™7	
Cat. Number/Size	T7-319-T025	25 tests
	T7-319-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	IgG (Fab)
Clone	4A11
Format	PE-Cy™7
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG2b
Specificity	The antibody 4A11 reacts with Fab fragment of human IgG.
Other names	Immunoglobulin G Fab fragment
Immunogen	Purified human IgG.

Preparation	The purified antibody is conjugated with tandem dye PE-Cy™7 under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

Images and References	www.exbio.cz
-----------------------	--

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.

Cy™ and CyDye™ are registered trademarks of GE Healthcare.