

Technical Data Sheet

Product	<b>Anti-HDAC6 Purified</b>	
Cat. Number/Size	<b>11-842-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	HDAC6
Clone	178
Format	Purified
Reactivity	Human
Application	ICC, ELISA, IP
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 178 recognizes zinc finger domain of human histone deacetylase 6 (HDAC6; an intracellular antigen). Crossreactivity with other species was not tested.
Other names	Histone deacetylase 6
Immunogen	human HDAC6
Entrez Gene ID	10013
Gene name	HDAC6
NCBI Full Gene Name	histone deacetylase 6
UniProt ID	Q9NSW6

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.