

Technical Data Sheet

Product	Anti-Hu IgM HRP	
Cat. Number/Size	1X-320-C025	0.025 mg
	1X-320-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	IgM
Clone	CH2
Format	HRP
Reactivity	Human
Application	WB, ELISA
Isotype	Mouse IgG1
Specificity	The antibody CH2 reacts with Fc fragment of human IgM.
Other names	immunoglobulin M
Immunogen	Purified human IgM.

Concentration	1 mg/ml
Preparation	The purified antibody is conjugated with horseradish peroxidase (HRP) of high specific activity and RZ=3.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing thimerosal (0.01%)
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.