



	Technical Data Sheet	
Product	Anti-Hu CD369 Purified	
Cat. Number/Size	11-874-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD369	
Clone	15E2	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IHC-F (reported)	
Application details	Flow cytometry: Recommended dilution: 2-5 µg/ml.	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody 15E2 recognizes an extracellular epitope LWEDGSTFSSN of human CD369 (Dectin-1), a 33 kDa transmembrane glycoprotein expressed predominantly on dendritic cells.	
Other names	Dectin-1, CLEC7A, BGR, CANDF4, SCARE2	
Workshop	HLDA X	
Immunogen	Extracellular part of human CD369 with hIgG FC tag	
Entrez Gene ID	64581	
Gene name	CLEC7A	
NCBI Full Gene Name	C-type lectin domain containing 7A	
UniProt ID	Q9BXN2	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.