

Technical Data Sheet

Product

**Anti-Kinesin Purified**

Cat. Number/Size

**11-434-C025**                      **0.025 mg**  
**11-434-C100**                      **0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	Kinesin
Clone	PAb (434)
Format	Purified
Reactivity	Mouse, Pig, Human
Application	WB, ICC
Isotype	Rabbit polyclonal
Specificity	The polyclonal antibody detects total level of endogenous kinesin protein (intracellular antigen).
Other names	KIN
Immunogen	stalk domain of human kinesin (aa 331-906) expressed in <i>E. coli</i> (FKHC3)

---

Concentration	1 mg/ml
Preparation	Purified by ligand affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References      [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.