



|                       |                                                                                                                                                                                                                                                                                                                                                                                                           |                  |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                       | Technical Data Sheet                                                                                                                                                                                                                                                                                                                                                                                      |                  |
| Product               | <b>Anti-Hu CD36 FITC</b>                                                                                                                                                                                                                                                                                                                                                                                  |                  |
| Cat. Number/Size      | <b>1F-451-T025</b>                                                                                                                                                                                                                                                                                                                                                                                        | <b>25 tests</b>  |
|                       | <b>1F-451-T100</b>                                                                                                                                                                                                                                                                                                                                                                                        | <b>100 tests</b> |
|                       | <b>For Research Use Only.</b>                                                                                                                                                                                                                                                                                                                                                                             |                  |
|                       | <b>Not for use in diagnostic or therapeutic procedures.</b>                                                                                                                                                                                                                                                                                                                                               |                  |
| Antigen               | CD36                                                                                                                                                                                                                                                                                                                                                                                                      |                  |
| Clone                 | TR9                                                                                                                                                                                                                                                                                                                                                                                                       |                  |
| Format                | FITC                                                                                                                                                                                                                                                                                                                                                                                                      |                  |
| Reactivity            | Human                                                                                                                                                                                                                                                                                                                                                                                                     |                  |
| Application           | FC (QC tested)                                                                                                                                                                                                                                                                                                                                                                                            |                  |
| Application details   | Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.                                                                                                                                                                                 |                  |
| Excitation laser      | blue (488 nm)                                                                                                                                                                                                                                                                                                                                                                                             |                  |
| Isotype               | Mouse IgG1                                                                                                                                                                                                                                                                                                                                                                                                |                  |
| Specificity           | The antibody TR9 reacts with an extracellular epitope of CD36 (GPIIb), a 85 kDa integral membrane glycoprotein expressed on platelets, macrophages, endothelial cells, early erythroid cells and megakaryocytes. The antibody TR9 cross-blocks binding of FITC-labeled standard antibody OKM5. Anti-CD36 antibodies inhibit adhesive functions (e.g. adherence of infected erythrocytes to target cells). |                  |
| Other names           | GPIIb, GPIV, PAS-4, FAT, Thrombospondin receptor, PASIV                                                                                                                                                                                                                                                                                                                                                   |                  |
| Immunogen             | Platelets                                                                                                                                                                                                                                                                                                                                                                                                 |                  |
| Entrez Gene ID        | 948                                                                                                                                                                                                                                                                                                                                                                                                       |                  |
| Gene name             | CD36                                                                                                                                                                                                                                                                                                                                                                                                      |                  |
| NCBI Full Gene Name   | CD36 molecule                                                                                                                                                                                                                                                                                                                                                                                             |                  |
| UniProt ID            | P16671                                                                                                                                                                                                                                                                                                                                                                                                    |                  |
| Preparation           | Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.                                                                                                                                                                                                             |                  |
| Formulation           | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                                                                                                                                                                                                                                                                                   |                  |
| Storage and handling  | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.                                                                                                                                                                                                                                                                                                                                  |                  |
| Images and References | <a href="http://www.exbio.cz">www.exbio.cz</a>                                                                                                                                                                                                                                                                                                                                                            |                  |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.