



	Technical Data Sheet	
Product	Anti-Por Transferrin Purified	
Cat. Number/Size	11-475-C025	0.025 mg
	11-475-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Transferrin	
Clone	PTF-02	
Format	Purified	
Reactivity	Pig	
Negative species	Horse, Sheep	
Application	ELISA, WB	
Isotype	Mouse IgG1	
Specificity	The antibody PTF-02 recognizes porcine serum transferrin, a 77 kDa single polypeptide chain glycoprotein (member of the iron binding family of proteins). It is synthesised in the liver and consists of two domains each having a high affinity reversible binding site for Fe3+.	
Other names	TFQTL1, PRO1557, PRO2086	
Immunogen	Porcine transferrin	
Entrez Gene ID	396996	
Gene name	TF	
NCBI Full Gene Name	transferrin	
UniProt ID	P09571	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.