

Technical Data Sheet

Product	Anti-Granzyme B Purified	
Cat. Number/Size	11-634-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Granzyme B
Clone	CLB-GB11
Format	Purified
Reactivity	Non-human primates, Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody CLB-GB11 recognizes granzyme B, a 31 kDa serine protease expressed intracellularly in activated Tc cells and NK cells.
Other names	GZMB, HLP, CTLA1, SECT
Immunogen	Human NK cell line YT-INDY-derived granzyme B
Entrez Gene ID	3002
Gene name	GZMB
NCBI Full Gene Name	granzyme B
UniProt ID	P10144

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.