

Technical Data Sheet

Product

Anti-Drebrin Purified

Cat. Number/Size

11-694-C025

0.025 mg

11-694-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	Drebrin
Clone	DBN-N-03
Format	Purified
Reactivity	Human
Application	FC (QC tested), WB, MC (CyTOF)
Application details	Western blotting: Positive control: REH, Nalm-6, HEK293 cells. Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody DBN-N-03 recognizes drebrin, an approximately 100-125 kDa intracellular regulator of actin cytoskeleton.
Other names	DBN, D0S117E
Immunogen	Bacterially expressed N-terminal fragment of recombinant human drebrin (aa 1-326)
Entrez Gene ID	1627
Gene name	DBN1
NCBI Full Gene Name	drebrin 1
UniProt ID	Q16643

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.