



	Technical Data Sheet	
Product	Anti-Ms CD11b Purified	
Cat. Number/Size	11-595-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD11b	
Clone	M1/70	
Format	Purified	
Reactivity	Mouse, Rabbit, Human, Non-human primates	
Application	FC (QC tested), IP, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1 µg/ml.	
Isotype	Rat IgG2b	
Specificity	The rat monoclonal antibody M1/70 detects an extracellular epitope of CD11b (integrin alphaM subunit), a type I transmembrane protein mainly expressed on monocytes/macrophages, granulocytes and NK-cells, which associates with CD18 to form Mac-1 integrin that plays important role in cell-cell interactions.	
Other names	Mac-1, Integrin alpha M, ITGAM, CR3A, MO1A, MAC1A	
Immunogen	B10 mouse spleen cells enriched for T cells	
Entrez Gene ID	16409	
Gene name	Itgam	
NCBI Full Gene Name	integrin subunit alpha M	
UniProt ID	G5E8F1	
Concentration	1 mg/ml	
Preparation	Purified by protein-G affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.