

Technical Data Sheet

Product	Anti-Hu CD18 Biotin	
Cat. Number/Size	1B-216-C025	0.025 mg
	1B-216-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD18
Clone	MEM-48
Format	Biotin
Reactivity	Human
Negative species	Dog
Application	FC (QC tested), IP, WB, IHC(P), IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody MEM-48 recognizes an extracellular epitope involving residues 534-546 in cysteine-rich repeat 3 of the CD18 antigen (integrin beta2 subunit; beta2 integrin). CD18 is a 90-95 kDa type I transmembrane protein expressed on all leukocytes.
Other names	Integrin beta2, Complement receptor C3 subunit beta, ITGB2
Immunogen	Leukocytes of a patient suffering from a LGL-type leukemia.
Entrez Gene ID	3689
Gene name	ITGB2
NCBI Full Gene Name	integrin subunit beta 2
UniProt ID	P05107

Concentration	1 mg/ml
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.