

Technical Data Sheet

Product	<b>Anti-LCK Purified</b>	
Cat. Number/Size	<b>11-269-C025</b>	<b>0.025 mg</b>
	<b>11-269-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	LCK
Clone	LCK-01
Format	Purified
Reactivity	Human
Negative species	Mouse
Application	FC, IP, WB, ICC
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Specificity	The antibody LCK-01 recognizes defined epitope (aa 22-36) of Lck, a 56 kDa Src-family protein tyrosine kinase (intracellular antigen).
Other names	p56Lck
Immunogen	Peptide corresponding to amino acids 22-36 in the sequence of human Lck.
Entrez Gene ID	3932
Gene name	LCK
NCBI Full Gene Name	LCK proto-oncogene, Src family tyrosine kinase
UniProt ID	P06239

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.