

Technical Data Sheet

Product	Anti-ABRA1 Purified	
Cat. Number/Size	11-619-C025	0.025 mg
	11-619-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	ABRA1
Clone	ABRA1-01
Format	Purified
Reactivity	Human
Application	FC, IP, WB
Isotype	Mouse IgG
Specificity	The mouse monoclonal antibody ABRA1-01 recognizes human ABRA1 (Abraxas, CCDC98), an adaptor protein involved in DNA repair, which migrates as a 45 kDa band on PAAGE under reducing conditions.
Other names	ABRA1, Abraxas 1, CCDC98, FAM175A
Immunogen	Recombinant protein corresponding to amino acids 1-313 of ABRA1 with N-terminal His6 tag
Entrez Gene ID	84142
Gene name	ABRAXAS1
NCBI Full Gene Name	abraxas 1, BRCA1 A complex subunit
UniProt ID	Q6UWZ7

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.